



THE  
SAN DIEGO  
CONFERENCE

## THE SAN DIEGO CONFERENCE THE GENETIC REVOLUTION

Christine L. Brakel, Ph. D., Senior Scientist

CORRESPONDENCE CONTACT: NAME AND TITLE

Enzo Diagnostics, Inc

516-694-7070

INSTITUTION

TELEPHONE

60 Executive Blvd.

516-694-7501

ADDRESS

FAX

Farmingdale, NY 11735

Hepatitis B DNA Detection

ADDRESS

POSTER TOPIC

FOR INTERNAL USE ONLY:

TEMP. #

PERM. #

POSTER ABSTRACT

COLORIMETRIC MICROPLATE HYBRIDIZATION ASSAYS FOR THE DETECTION OF HEPATITIS B VIRUS DNA IN SERUM. Christine L. Brakel, Enzo Diagnostics, Inc. An antibody-enhanced microplate hybridization assay (EMHA) was used for the direct detection of Hepatitis B Virus (HBV) DNA in serum. Eighteen of 60 (30%) Hepatitis B Surface Antigen (HBsAg) positive specimens contained HBV DNA while only 2 of 372 (0.5%) HBsAg-negative specimens were positive for HBV DNA. Forty-six specimens (37 HBsAg positive and 9 HBsAg negative) and 4 titration standards for the EMHA were tested independently for HBV DNA using the Abbott radiological Hepatitis B Viral DNA Hybridization test. Of 9 HBsAg-positive specimens that were positive for HBV DNA by the EMHA, 7 were positive by the Abbott test. The Abbott test identified 2 other HBV DNA-positive specimens (>3pg/ml DNA) as well as 2 specimens that contained equivocal amounts of HBV DNA (2 and 3.4pg/ml). Correlating the results of these two quantitative DNA hybridization tests indicated a detection limit for the EMHA of 15 to 20 pg of HBV DNA.

In a separate study, the EMHA was again compared to the Abbott HBV DNA test and to the results of polymerase chain reaction (PCR) amplification of DNA isolated from serum. Those results confirmed a close correlation between the results of the EMHA and the Abbott test. As expected, both direct detection methods failed to detect several samples that contained HBV DNA when the nucleic acid was subjected to PCR amplification prior to assay in a standard Microplate Hybridization Assay. Unexpectedly, the results of PCR amplification indicated that in some specimens, there had been gene mutation and rearrangement in the Core antigen sequences of the Hepatitis B genome.

1988 ASM ANNUAL MEETING  
Miami Beach, Fla. 8-13 May 1988

Poster

Official Abstract Form

(Read all instructions before typing)

Direct Detection of HSV Infection by DNA Hybridization,  
BRAKEL, C.L., POLLICE, M.A., and VICTOR, J. ENZO  
Biochem Inc., New York, N.Y.

The *ColorGene*<sup>TM</sup> DNA Hybridization Test for HSV Confirmation was applied to the direct detection of HSV-infected cells from herpes lesions. For this clinical trial 2 swabs from suspected herpes lesions were taken. The first swab specimen was used for culture isolation of HSV and material from the second swab was smeared onto glass microscope slides, fixed, and analyzed by *in situ* DNA hybridization. The results of this trial, carried out at 4 clinical virology laboratories, were 184 culture +/DNA +; 102 culture -/DNA -; 30 culture +/DNA -; and 24 culture -/DNA + (sensitivity = 86%, specificity = 81%). Of the 24 false positive specimens, 14 were from patients with a history of herpes infection, suggesting that the virus in these specimens had not grown in culture. On the average, specimens that tested as false negatives took longer to develop cytopathic effect in culture than specimens that tested true positive. The sensitivity of the test for vesicles, pustules, ulcers, and crusts was 92%, 87%, 82%, and 71%, respectively, indicating that direct testing is optimal for early as opposed to healing lesions, as is also true for isolation of virus in culture. Although the test does not distinguish HSV 1 from HSV 2, the test was equally sensitive for direct detection of both HSV types. The test is suitable for use in clinics and/or physicians' offices as the method requires no sophisticated equipment, is easy to interpret, and results can be observed with a simple light microscope. After fixation of the specimen (5-10 minutes), the test requires only 45-60 minutes until results can be read.

Instructions

Indicate below the subject category designation from the list on p. iv, check your poster or slide session preference, complete the check list on the reverse side of this sheet, and sign your name in the space provided.

Indicate category designation from page iv.

S-4

Category designation

Poster/Slide Session Preference

Because of the flexibility in programming afforded by poster sessions, the Program Committee will attempt to schedule *all* abstracts which (i) are considered by elected divisional officers to be of acceptable quality and (ii) conform to rules established by the Program Committee. The decision on whether an abstract is scheduled in a slide or a poster session will be made by the elected Program Committee, which will be guided (but not bound) by the preference of the authors. Approximately 75% of the abstracts will be scheduled in poster sessions. By submitting an abstract, the author agrees that the paper will be presented as scheduled.

Please check one: ☒ Poster session preferred ☐ Slide session preferred ☐ No preference

Please provide telephone number of signing author 212) 741-3838 x 115

Area code

Tues AM

#166