

Slide Session

Hard Copy  
of Slides  
Attached

1988 ASM ANNUAL MEETING  
Miami Beach, Fla. 8-13 May 1988

Official Abstract Form

(Read all instructions before typing)

An Automatable, Colorimetric DNA Hybridization Test for *M. tuberculosis* Confirmation, BRAKEL, C.L., DONEGAN, J.J., LINN, C-I.P., MOLINA, M., POLLICE, M.A., WANG, Z., and YANG, H. L. ENZO Biochem Inc., New York, N. Y.

An oligonucleotide-based DNA hybridization test for confirmation of *M. tuberculosis* (MTB) cultures has been developed that is amenable to either partial or complete automation. Following lysis (10 min.) of cultured specimens, the hybridization is carried out in two steps and can be accomplished in less than 2 hours (20-30 minutes "hands on" time), even when as many as 30-60 specimens are to be analyzed. The lysed cultured specimens are first hybridized against one modified oligomeric probe in solution and are then allowed to hybridize to a second probe coated onto wells of microtiter (ELISA) plates. After hybridization and washing, the hybrids are detected with streptavidin-biotinylated horseradish peroxidase. Signal is generated by enzymatic conversion of hydrogen peroxide and o-phenylenediamine. Results can be read by eye, or quantitated with an ordinary ELISA photometer. To date the test has been 100% sensitive and specific. In a blind confirmation of 86 clinical isolates, 64 were correctly identified as MTB and 22 as non-MTB. In addition, 83 different species of bacteria, including 22 species of mycobacteria have been identified correctly as non-MTB. This methodology is suitable for the automated confirmation of any cultured organism provided suitable probes are available.

Instructions

Indicate below the subject category designation from the list on p. iv, check your poster or slide session preference, complete the check list on the reverse side of this sheet, and sign your name in the space provided.

Indicate category designation from page iv.

U-1

Category designation

Poster/Slide Session Preference

Because of the flexibility in programming afforded by poster sessions, the Program Committee will attempt to schedule *all abstracts* which (i) are considered by elected divisional officers to be of acceptable quality and (ii) conform to rules established by the Program Committee. The decision on whether an abstract is scheduled in a slide or a poster session will be made by the elected Program Committee, which will be guided (but not bound) by the preference of the authors. Approximately 75% of the abstracts will be scheduled in poster sessions. By submitting an abstract, the author agrees that the paper will be presented as scheduled.

Please check one: ☐ Poster session preferred ☒ Slide session preferred ☐ No preference

Please provide telephone number of signing author (212) 741-3838 x 125

Area code

Slide Session 256(u) Thurs AM  
11:20 11:20

**An Automatable, Colorimetric  
DNA Hybridization Test for  
M. tuberculosis Confirmation**

**Special Thanks to**

**Jim Donegan**

**Patsy Lin**

**Margarita Molina**


**Marjorie Pollice**

**Zwang Wang**

**Huey Lang Yang**


## Map and Sequences from MTB probe p24861

p24861 is a 7 kb insert in the Bam H1 site of pIB1 76, grown in HB101




Oligonucleotides 4, 5, 7 and 8 are 30 base sequences synthesized using an Applied Biosystems Synthesizer. The two pairs are from opposite strands of the DNA.

# Specimen Lysis and First Hybridization




## Second Hybridization and Wash



# Detection of Hybridized MTB DNA

7



(NOTE: Sequences are examples only, and are not the actual sequences used.)

## Cross Reaction Studies

Acinetobacter calcoaceticus	Enterobacter aerogenes	R. sputi
A. Iwoffii	Escherichia coli	Rhodospirillum rubrum
Actinomadura madurae	Fusobacterium nucleatum	Staphylococcus aurea
Actinoplanes italicus	Haemophilus influenzae	Streptococcus mitis
Arthrobacter oxydans	Klebsiella pneumoniae	S. pneumoniae
Bacillus subtilis	Legionella pneumophila	Vibrio parahaemolyticus
Bacterionema matruchotii	Microbacterium lacticum	Yersinia enterocolitica
Bacteroides fragilis	Mycoplasma hominis	
Branhamella Catarrhalis	M. pneumoniae	
Brevibacterium linens	Neisseria gonorrhea	
Campylobacter jejuni	N. lactamica	
Chromobacterium violaceum	N. meningitidis	
Clostridium perfringens	Nocardia asteroides	
Corynebacterium aquaticum	N. brasiliensis	
C. diphtheriae	N. otitis-caviarum	
C. genitalium	Nocardiosis dassonvillei	
C. haemolyticum	Oerskovia turbata	
C. minutissimum	O. xanthineolytica	
C. pseudodiphtheriticum	Propionibacterium acnes	
C. pseudogenitalium	Proteus mirabilis	
C. pseudotuberculosis	Pseudomonas aeruginosa	
C. pyogenes	P. cepacia	
C. renale	Rahnella aquatilis	
C. striatum	Rhodococcus aichiensis	
C. xerosis	R. aurantiacus	
Deinococcus radiodurans	R. bronchialis	
Dermatophilus congolensis	R. chubeuensis	
Derxia gummosa	R. equi	

**63 Bacterial  
Species found  
to be Not  
Cross Reactive**

## Cross Reaction Studies Mycobacterial Species

### POSITIVE REACTIONS

#### MTB Complex

*Mycobacterium africanum*

*Mycobacterium bovis*

*Mycobacterium bovis* BCG

*Mycobacterium tuberculosis*

*Mycobacterium microti*

Positive for MTB Complex  
Negative for 25 other  
Mycobacterial Species

### NEGATIVE REACTIONS

*M. asiaticum*

*M. avium*

*M. chelonae*

*M. flavescens*

*M. fortuitum*

*M. gastri*

*M. gordonae*

*M. haemophilum*

*M. intracellulare*

*M. kansasii*

*M. malmoense*

*M. marinum*

*M. nonchromogenicum*

*M. phlei*

*M. scrofulaceum*

*M. shimoidei*

*M. simiae*

*M. smegmatis*

*M. szulgai*

*M. terrae*

*M. thermoresistibile*

*M. triviale*

*M. ulcerans*

*M. vaccae*

*M. xenopi*



## Results of Testing

Study 1

Study 2

Culture

Culture/GenProbe

DNA Assay	Culture	
	+	-
+	72	1 *
-	1 * *	49

DNA Assay	Culture/GenProbe	
	+	-
+	55	0
-	0	38

+

+

DNA Assay

DNA Assay

-

-

\*Originally identified as *M. xenopi*.  
Later found to contain MTB and  
*Corynebacterium pseudotuberculosis*.

All specimens were originally  
identified by use of the GenProbe  
MTB complex confirmation test.

\*\*Originally identified as MTB.  
Later found to be *Brevibacterium*  
*linens*.