

Slide Session

Hard Copy
of Slides
Attached

1988 ASM ANNUAL MEETING
Miami Beach, Fla. 8-13 May 1988

Official Abstract Form

(Read all instructions before typing)

An Automatable, Colorimetric DNA Hybridization Test for *M. tuberculosis* Confirmation, BRAKEL, C.L., DONEGAN, J.J., LINN, C-I.P., MOLINA, M., POLLICE, M.A., WANG, Z., and YANG, H. L. ENZO Biochem Inc., New York, N. Y.

An oligonucleotide-based DNA hybridization test for confirmation of *M. tuberculosis* (MTB) cultures has been developed that is amenable to either partial or complete automation. Following lysis (10 min.) of cultured specimens, the hybridization is carried out in two steps and can be accomplished in less than 2 hours (20-30 minutes "hands on" time), even when as many as 30-60 specimens are to be analyzed. The lysed cultured specimens are first hybridized against one modified oligomeric probe in solution and are then allowed to hybridize to a second probe coated onto wells of microtiter (ELISA) plates. After hybridization and washing, the hybrids are detected with streptavidin-biotinylated horseradish peroxidase. Signal is generated by enzymatic conversion of hydrogen peroxide and o-phenylenediamine. Results can be read by eye, or quantitated with an ordinary ELISA photometer. To date the test has been 100% sensitive and specific. In a blind confirmation of 86 clinical isolates, 64 were correctly identified as MTB and 22 as non-MTB. In addition, 83 different species of bacteria, including 22 species of mycobacteria have been identified correctly as non-MTB. This methodology is suitable for the automated confirmation of any cultured organism provided suitable probes are available.

Instructions

Indicate below the subject category designation from the list on p. iv, check your poster or slide session preference, complete the check list on the reverse side of this sheet, and sign your name in the space provided.

Indicate category designation from page iv.

U-1

Category designation

Poster/Slide Session Preference

Because of the flexibility in programming afforded by poster sessions, the Program Committee will attempt to schedule *all abstracts* which (i) are considered by elected divisional officers to be of acceptable quality and (ii) conform to rules established by the Program Committee. The decision on whether an abstract is scheduled in a slide or a poster session will be made by the elected Program Committee, which will be guided (but not bound) by the preference of the authors. Approximately 75% of the abstracts will be scheduled in poster sessions. By submitting an abstract, the author agrees that the paper will be presented as scheduled.

Please check one: ☐ Poster session preferred ☒ Slide session preferred ☐ No preference

Please provide telephone number of signing author (212) 741-3838 x 125

Area code

Slide Session 256(u) Thurs AM
11:20 11:25

**An Automatable, Colorimetric
DNA Hybridization Test for
M. tuberculosis Confirmation**

Special Thanks to

Jim Donegan

Patsy Lin

Margarita Molina


Marjorie Pollice

Zwang Wang

Huey Lang Yang


Map and Sequences from MTB probe p24861

p24861 is a 7 kb insert in the Bam H1 site of pIB1 76, grown in HB101




Oligonucleotides 4, 5, 7 and 8 are 30 base sequences synthesized using an Applied Biosystems Synthesizer. The two pairs are from opposite strands of the DNA.

Specimen Lysis and First Hybridization




Second Hybridization and Wash



Detection of Hybridized MTB DNA

7



(NOTE: Sequences are examples only, and are not the actual sequences used.)

Cross Reaction Studies

| | | |
|-----------------------------|--------------------------|-------------------------|
| Acinetobacter calcoaceticus | Enterobacter aerogenes | R. sputi |
| A. Iwoffii | Escherichia coli | Rhodospirillum rubrum |
| Actinomadura madurae | Fusobacterium nucleatum | Staphylococcus aurea |
| Actinoplanes italicus | Haemophilus influenzae | Streptococcus mitis |
| Arthrobacter oxydans | Klebsiella pneumoniae | S. pneumoniae |
| Bacillus subtilis | Legionella pneumophila | Vibrio parahaemolyticus |
| Bacterionema matruchotii | Microbacterium lacticum | Yersinia enterocolitica |
| Bacteroides fragilis | Mycoplasma hominis | |
| Branhamella Catarrhalis | M. pneumoniae | |
| Brevibacterium linens | Neisseria gonorrhea | |
| Campylobacter jejuni | N. lactamica | |
| Chromobacterium violaceum | N. meningitidis | |
| Clostridium perfringens | Nocardia asteroides | |
| Corynebacterium aquaticum | N. brasiliensis | |
| C. diphtheriae | N. otitis-caviarum | |
| C. genitalium | Nocardiosis dassonvillei | |
| C. haemolyticum | Oerskovia turbata | |
| C. minutissimum | O. xanthineolytica | |
| C. pseudodiphtheriticum | Propionibacterium acnes | |
| C. pseudogenitalium | Proteus mirabilis | |
| C. pseudotuberculosis | Pseudomonas aeruginosa | |
| C. pyogenes | P. cepacia | |
| C. renale | Rahnella aquatilis | |
| C. striatum | Rhodococcus aichiensis | |
| C. xerosis | R. aurantiacus | |
| Deinococcus radiodurans | R. bronchialis | |
| Dermatophilus congolensis | R. chubeuensis | |
| Derxia gummosa | R. equi | |

**63 Bacterial
Species found
to be Not
Cross Reactive**

Cross Reaction Studies Mycobacterial Species

POSITIVE REACTIONS

MTB Complex

Mycobacterium africanum
Mycobacterium bovis
Mycobacterium bovis BCG
Mycobacterium tuberculosis
Mycobacterium microti

Positive for MTB Complex
Negative for 25 other
Mycobacterial Species

NEGATIVE REACTIONS

| | | |
|--------------------------|----------------------------|-----------------------------|
| <i>M. asiaticum</i> | <i>M. kansasii</i> | <i>M. szulgai</i> |
| <i>M. avium</i> | <i>M. malmoense</i> | <i>M. terrae</i> |
| <i>M. chelonae</i> | <i>M. marinum</i> | <i>M. thermoresistibile</i> |
| <i>M. flavescens</i> | <i>M. nonchromogenicum</i> | <i>M. triviale</i> |
| <i>M. fortuitum</i> | <i>M. phlei</i> | <i>M. ulcerans</i> |
| <i>M. gastri</i> | <i>M. scrofulaceum</i> | <i>M. vaccae</i> |
| <i>M. gordonae</i> | <i>M. shimoidei</i> | <i>M. xenopi</i> |
| <i>M. haemophilum</i> | <i>M. simiae</i> | |
| <i>M. intracellulare</i> | <i>M. smegmatis</i> | |

Results of Testing

Study 1

Study 2

Culture

Culture/GenProbe

| DNA Assay | Culture | |
|-----------|---------|-----|
| | + | - |
| + | 72 | 1 * |
| - | 1 * * | 49 |

| DNA Assay | Culture/GenProbe | |
|-----------|------------------|----|
| | + | - |
| + | 55 | 0 |
| - | 0 | 38 |

*Originally identified as *M. xenopi*.
Later found to contain MTB and
Corynebacterium pseudotuberculosis.

All specimens were originally
identified by use of the GenProbe
MTB complex confirmation test.

**Originally identified as MTB.
Later found to be *Brevibacterium*
linens.